



# ARX

## 4.0 Cannulated Screw System



## Contents

3	Indications
4	Features and Benefits
5	Ordering Information
	- 4.0 Cannulated Screws
6	Ordering Information
	- Instrument
7	Set Configuration Suggestions

## Indications

### **Fixation of Small Bone Fragments**

- Fracture of the phalanges, metacarpals, carpals and distal radius
- Fractures of the phalanges, metatarsals, tarsals and ankle

### **Arthrodeses of Corrective Osteotomies of Small Bones**

- Arthrodeses of small joints of hand and foot
- Bunionectomies and Osteotomies

## Thread Length Choice

For best fit of distal threads into far bone fragment to bear increased interfragmentary compression

## Flutes

- Self-Tapping Flutes for Reduced OR time by eliminating the need for tapping
- Reverse Cutting Flute for ease of screw removal

## Shaft and Core diameter equal for increased strength of screw

## Low Profile Screw Head

Reduced Possibility of Soft Tissue Irritation

## Titanium Alloy Material

- Increased Compatibility with CT/MRI
- Biocompatible, Corrosion-resistant and Non-toxic in a Biological Environment

## Screw Options


Thread Length(Th) : 7mm ~ 23mm

Screw Length(L) : 20mm ~ 70mm

Outer Diameter : 4.0mm



Cannulated Screw(Starix)



ARIX Cannulated Screws

[Unit: mm]

No.	Code	Head Diameter(A)	Thread Diameter(B)	Thread Length(Th)	Screw Length(L)
1	40-CO-020-T	Ø5.0	Ø4.0	7	20
2	40-CO-022-T			7	22
3	40-CO-024-T			8	24
4	40-CO-026-T			9	26
5	40-CO-028-T			9	28
6	40-CO-030-T			10	30
7	40-CO-034-T			11	34
8	40-CO-038-T			12	38
9	40-CO-042-T			14	42
10	40-CO-046-T			15	46
11	40-CO-050-T			16	50
12	40-CO-054-T			18	54
13	40-CO-058-T			19	58
14	40-CO-062-T			21	62
15	40-CO-066-T			22	66
16	40-CO-070-T			23	70



## 112-40-701

Drill Bit for Ø4.0mm Cannulated Screw(AO Type)



## 113-HF-617

T15 Driver Inner Hole Ø1.4mm



## 111-072

Cannulated Countersink Drill Bit



## 111-063

Screwdriver Handle



## 111-062

Depth Gauge for Cannulated Screw




## 111-098

Drill Guide for Ø4.0mm Cannulated Screw



## 111-068-5

Guide Pin Ø1.3mm



## 111-067-5

Stylet Ø1.3mm



## 111-096

Dispenser

## Set Configuration



### 112-120

ARIX 4.0 Cannulated Screw System Tray

NO	Code	Description	Set Quantity
<b>Screws</b>			
1	40-CO-020-T	Length 20mm	4
2	40-CO-022-T	Length 22mm	4
3	40-CO-024-T	Length 24mm	4
4	40-CO-026-T	Length 26mm	4
5	40-CO-028-T	Length 28mm	4
6	40-CO-030-T	Length 30mm	4
7	40-CO-034-T	Length 34mm	4
8	40-CO-038-T	Length 38mm	4
9	40-CO-042-T	Length 42mm	4
10	40-CO-046-T	Length 46mm	4
11	40-CO-050-T	Length 50mm	4
12	40-CO-054-T	Length 54mm	4
13	40-CO-058-T	Length 58mm	4
14	40-CO-062-T	Length 62mm	4
15	40-CO-066-T	Length 66mm	4
16	40-CO-070-T	Length 70mm	4

NO	Code	Description	Set Quantity
<b>Instruments</b>			
1	112-120	ARIX 4.0 Cannulated Screw System Tray	1
2	112-40-701	Drill Bit for Ø4.0mm Cannulated Screw (AO Type)	2
3	113-HF-617	T15 Driver Inner Hole Ø1.4mm	2
4	111-072	Cannulated Countersink Drill Bit	1
5	111-063	Screwdriver Handle	1
6	111-062	Depth Gauge for Cannulated Screw	1
7	111-098	Drill Guide for Ø4.0mm Cannulated Screw	1
8	111-068-5	Guide Pin Ø1.3mm	10
9	111-067-5	Stylet Ø1.3mm	2
10	111-096	Dispenser	1



## JEIL MEDICAL CORPORATION

702•703•704•705•706•804•805•807•812-ho, 55, Digital-ro 34-gil, Guro-gu, Seoul, 08378, Rep. of Korea

**INTERNATIONAL SALES\_** T 82.2.850.3935 F 82.2.850.3537 E [global@jeilmed.co.kr](mailto:global@jeilmed.co.kr)

**DOMESTIC SALES\_** T 82.2.850.3500 F 82.2.850.3535 E [jeil@jeilmed.co.kr](mailto:jeil@jeilmed.co.kr)